

SPECIFICATION

SYSTEMS AND METHODS TO IDENTIFY AND DISABLE RE-USED SINGLE USE DEVICES BASED ON DETECTING ENVIRONMENTAL CHANGES

FIELD OF THE INVENTION

5 The present invention is directed to medical systems, and more specifically, to systems, apparatus and methods for limiting the re-usage of medical probes, such as catheters and surgical probes.

BACKGROUND

10 Catheters, surgical probes and related probe devices (collectively referred to, hereinafter, as "medical probes" or "probes") are used today in diagnostic and therapeutic medical procedures that require surgical or minimally invasive access to target tissue areas within interior regions of the body. During these procedures, a physician locates the distal end of the medical probe at the target site by, in the case of a catheter, steering the medical
15 probe through a main vein or artery (typically, the femoral vein or artery), or, in the case of a surgical probe, advancing the medical probe through the surgical opening leading to the target site.

 The physician then operates the medical probe to activate a functional component mounted at the probe distal end, thereby providing therapeutic treatment to and/or
20 diagnosis of the interior region. Due to the potential of passing any of a variety of dangerous diseases from one patient to another, prudent consideration dictates that the re-use of such probes should be prevented or, at the least, minimized. Typically, medical probes that are re-used are sterilized between uses to kill any disease-causing agents and to

remove any tissue that has collected on the medical probe during the previous use.

Sterilization of used medical probes, however, is not fool-proof, and oftentimes ineffective when tissue located on the medical probe is not exposed to the full effect of the sterilization process. Thus, even a sterilized medical probe may pose a threat to patients.

5 Minimizing re-use of medical probes that provide therapy becomes even more critical. During the therapeutic process, it is important for the physician to control the therapeutic component at the probe distal end carefully and precisely, so that adverse damage to a therapeutic component of the medical probe does not occur. Sterilization of and re-use of therapeutic medical probes subjects the therapeutic component to
10 mechanical, chemical and/or thermal stress, thus jeopardizing the control that the physician may have of the therapeutic component.

The need for careful and precise control over a therapeutic medical probe is especially critical during procedures that ablate tissue within the heart. These procedures, called electrophysiological therapy, are becoming more widespread for treating cardiac
15 rhythm disturbances. When inside the desired chamber of the heart, the physician manipulates a steering mechanism to place one or more electrodes located at the distal end of the medical probe into direct contact with the heart tissue to be ablated. The physician then directs radio frequency energy from the electrodes through the tissue to an indifferent electrode, thereby ablating the tissue and forming a lesion. If the electrodes or
20 the electrical wires connected thereto are worn or faulty, however, ablation may be ineffective and, in the worst case, may cause charring of the heart tissue.

Preventing or limiting re-usage of medical probes, while still allowing legitimate use of these probes, is made difficult by a possible scenario wherein the physician uses the

medical probe, temporarily disconnects the probe from the control unit, and reconnects the probe to the control unit to continue the procedure. Thus, there is a danger of deeming the continued procedure to be re-usage of the medical probe, which may result in the probe being rendered prematurely inoperable.

5 Thus, it would be desirable to provide an improved medical system for minimizing the re-usage of medical probes, while still allowing legitimate use of these probes.

SUMMARY OF THE INVENTION

The present invention is directed to improved apparatus and methods for limiting
10 the usage of medical probes, such as, e.g., catheters and surgical probes.

In a preferred method performed in accordance with a first aspect of the invention, the sterilization of a medical probe is monitored by detecting an environmental condition to which the medical probe is exposed, electronically storing a probe sterilization indicator with the medical probe if the detected environmental condition indicates exposure of the
15 medical probe to a sterilization cycle and determining the presence of the probe sterilization indicator. By way of non-limiting example, either temperature, moisture, pressure and/or a chemical is the detected environmental condition. The medical probe is conditionally operated based on the presence of the probe sterilization indicator. The operation of the medical probe may also be based on the presence of a probe usage
20 indicator, which indicates if the medical probe has been previously operated.

In another preferred method performed in accordance with a further aspect of the invention, the usage of a medical probe is limited by detecting an environmental condition to which the medical probe is exposed, electronically storing a probe sterilization indicator

in the medical probe if the detected environmental condition indicates exposure of the medical probe to a sterilization cycle, determining a presence of the medical probe indicator, and conditionally operating the medical probe based on the presence of the probe sterilization indicator. A probe usage indicator indicating that the medical probe has
5 been operated can be stored in the medical probe. Conditional operation of the medical probe can also be based on the probe usage indicator. Alternatively, an estimated probe sterilization date or probe manufacture date can be stored in the medical probe when the medical probe is manufactured, and an actual probe sterilization date can be stored in the medical probe upon actual sterilization of the medical probe. Conditional operation of the
10 medical probe can be based on the estimated probe sterilization date and the actual probe sterilization date. By way of non-limiting example, the medical probe can be prevented from being operated if the actual probe sterilization date is later than the estimated probe sterilization date. Optionally, the actual probe sterilization date can be used to enforce the shelf-life of the medical probe. By way of non-limiting example, operation of the medical
15 probe can be prevented if the difference between the actual probe sterilization date and a reference date, such as, e.g., the current date, is greater than a predetermined period of time.

In a preferred embodiment constructed in accordance with still another aspect of the invention, a control unit for connection to a medical probe having electronic storage
20 componentry, comprises control circuitry configured to electrically couple to the electronic storage componentry for reading data from the electronic storage componentry, and for conditionally operating the medical probe based on a presence of a probe sterilization indicator in the data. By way of non-limiting example, the preferred embodiment of the

control unit is an RF generator having a microprocessor as the control circuitry. The microprocessor is configured for conditionally operating the medical probe based on the presence of the probe sterilization indicator. The microprocessor can also be configured for conditionally operating the medical probe based on the presence of a probe usage indicator indicating previous operation of the medical probe. When the medical probe is operated, and more preferably, when the medical probe is effectively operated, the probe usage indicator is stored in the electronic storage componentry and the probe sterilization indicator is cleared.

In another preferred embodiment constructed in accordance with still a further aspect of the invention, a medical probe includes an elongate member carrying electronic storage componentry configured for detecting an environmental condition to which the medical probe is exposed and storing a probe sterilization indicator in the medical probe if the detected environmental condition indicates exposure of the medical probe to a sterilization cycle. By way of non-limiting example, in the preferred embodiment, the electronic storage component is a non-volatile storage component, such as, e.g., non-volatile RAM and a battery. Alternatively, the non-volatile storage components can comprise other components that have environmental sensing capabilities. For example, the non-volatile component can comprise an electrolytic sensor or chemfet.

Other and further objects, features, aspects, and advantages of the present invention will become better understood with the following detailed description of the accompanying drawings.

BRIEF DESCRIPTION OF THE DRAWINGS

The drawings illustrate both the design and utility of preferred embodiments of the present invention, in which:

Fig. 1 is a perspective view of an embodiment of a catheter system in accordance with various aspects of the present invention;

Fig. 2 is a schematic of a preferred embodiment of a medical probe system;

Fig. 3 is a flow diagram of a preferred method of conditionally operating a medical probe employed in the medical probe system of Fig. 2;

Fig. 4 is a schematic of an alternative preferred embodiment of a medical probe system;

Fig. 5 is a flow diagram of a preferred method of storing a probe sterilization indicator within a medical probe employed in the medical probe system of Fig. 4;

Fig. 6 is a flow diagram of a preferred method of conditionally operating a medical probe employed in the medical probe system of Fig. 4;

Fig. 7 is a schematic of another alternative preferred embodiment of a medical probe system;

Fig. 8 is a schematic of a catalog employed in the medical probe system of Fig. 7;

Fig. 9 is a flow diagram of a preferred method of obtaining historical operational information and conditionally operating the medical probe employed in the medical probe system of Fig. 7;

Fig. 10 is a schematic of still another alternative preferred embodiment of a medical probe system;

Fig. 11 is a flow diagram of a preferred method of conditionally operating a medical probe employed in the medical probe system of Fig. 10; and

Fig. 12 is a perspective view of an embodiment of a surgical probe system in accordance with various aspects of the present invention.

5

DETAILED DESCRIPTION OF THE PREFERRED EMBODIMENTS

Catheter System

Fig. 1 depicts a preferred embodiment of a catheter system 10 constructed in accordance with the present invention. The system 10 generally includes a radiofrequency generator 12, which delivers radio frequency energy, and a steerable catheter 14 carrying a radiofrequency emitting tip ablation electrode 16.

In the illustrated embodiment, the system 10 operates in a monopolar mode. In this arrangement, the system 10 includes a skin patch electrode that serves as an indifferent second electrode 18. In use, the indifferent electrode 18 is attached to the patient's back or other exterior skin area. Alternatively, the system 10 can be operated in a bipolar mode.

In this mode, the catheter 14 carries both electrodes. In the illustrated embodiment, the ablation electrode 16 and the indifferent electrode 18 are made of platinum. The catheter 14 carries a temperature sensor 19 (shown in Fig. 2) adjacent the ablation electrode 16 to measure the temperature to which the tissue is exposed to during ablation.

The catheter 14 further includes a handle 20, a guide tube 22, and a tip 24, which carries the ablation electrode 16. The handle 20 encloses a steering mechanism 26 for the catheter tip 24. A cable 28 extending from the rear of the handle 20 has plugs (not shown). The plugs connect the catheter 14 to the generator 12 for conveying

radiofrequency energy to the ablation electrode 16. The radiofrequency energy heats the tissue to form the lesion.

Left and right steering wires (not shown) extend through the guide tube 22 to interconnect the steering mechanism 26 to the left and right sides of the tip 24. Rotating the steering mechanism 26 to the left pulls on the left steering wire, causing the tip 24 to bend to the left. Also, rotating the steering mechanism 26 to the right pulls on the right steering wire, causing the tip 24 to bend to the right. In this way, the physician steers the ablation electrode 16 into contact with the tissue to be ablated.

Upon initial power on, the generator 12 is placed in a standby mode, which allows the physician to adjust the setpoint parameters. These setpoint parameters include the magnitude of the RF power delivered, tissue temperature and duration of RF power delivery.

The RF power delivered by the generator 12 is set using a power control input 30. The actual RF power delivered by the generator 12 is reported by a power meter 32.

During RF energy delivery, the generator 12 adjusts power output to maintain an actual measured temperature at the temperature setpoint. Thus, the measured power may be lower than the setpoint power depending on the measured tissue temperature and/or tissue impedance.

The desired temperature to which the ablated tissue is exposed is set using a temperature control input 34. The actual temperature to which the ablated tissue is exposed is reported by a temperature gauge 36. The desired duration of RF power application is set using a timer 38. A counter 40 tracks the elapsed time from initial delivery of RF power to the tissue, and counts from zero to the setpoint duration. When

loss of contact with tissue is detected, the counter 40 stops. Contact between the ablation electrode 16 and the tissue is measured with an impedance meter 42.

The generator 12 includes an RF power control button 44, the depression of which places the generator 12 in deliver mode. While in the deliver mode, the generator 12 delivers RF energy to the tissue in contact with the ablation electrode 16 until the count displayed by the counter 40 reaches the setpoint duration or until the RF power control button 44 is depressed again.

Medical Probe System

Fig. 12 depicts a preferred embodiment of a surgical probe system 200 constructed in accordance with the present invention. The surgical probe system 200 is similar to the catheter system 10 described above. Accordingly, elements of the surgical probe system 200 that are similar to those of the catheter system 10 have been given common reference numbers. Like the catheter system 10, the system 200 includes a radiofrequency generator 12, the details of which are described above. The system 200 further includes a surgical probe 202 that includes a relatively short, relatively stiff shaft 204, a handle 206, and a distal section 208. The shaft 204 may be from about 4 inches to 18 inches in length and is preferably about 6 to 8 inches. The distal section 208 may be from about 1 inch to 10 inches in length and is preferably about 4 to 6 inches. The shaft 204 supports a radiofrequency emitting tip ablation electrode 16, much like that described above, and a plurality of coil electrodes 210. This embodiment is particularly useful because it can be easily inserted into the patient through an introducing port such as a trocar.

Additional information concerning the above-described and other surgical probes may be found in U.S. Application Ser. No. 09/072,872, filed May 5, 1998, entitled "Surgical Method and Apparatus for Positioning a Diagnostic or Therapeutic Element Within the Body," which is incorporated herein by reference.

5 Interaction between the surgical probe 202 and the RF generator 12 is similar to that described above between the catheter 14 and the RF generator 12. For purposes of brevity, such operation will not be repeated.

Referring to Fig. 2, a medical probe system 300, which can take the form of the catheter system 10 or surgical probe system 200 described above, is depicted. A medical
10 probe 302 (i.e., the catheter 14 or surgical probe 202) includes a probe identification component 50, which is capable of storing a probe identification code C_i . The probe identification code C_i uniquely identifies the medical probe 302, such as by a serial number, or uniquely identifies a characteristic of the medical probe 302, such as, e.g., a particular performance and/or physical characteristic of the medical probe 302. In the
15 illustrated embodiment, the probe identification component 50 is a non-volatile storage component, and particularly is a non-solid state device, such as, e.g., one or more resistors, the resistance value of which represents the probe identification code C_i . Alternatively, the probe identification component 50 is a microchip, such as, e.g., a ROM, which has been pre-programmed with a digital identification code C_i .

20 The medical probe 302 further includes a probe usage limitation component 51. In this embodiment, the probe usage limitation component 51 is a non-volatile electronic storage component that is capable of storing data. Of course, the entire data need not be stored in a single component, but may be stored in several components without straying

from the principles taught by this invention. For purposes of illustration, the probe identification component 50 and the probe usage limitation component 51 are depicted as separate components. It should be noted, however, that the functions performed by the probe identification component 50 and the probe usage limitation component 51 can, in most instances, be performed within a single component.

Limiting Re-Use of Devices Based on Time Elapsed From First Use

The medical probe system 300 is capable of identifying and disabling re-used single use devices based on the time elapsed from the first therapeutic use of the medical probe 302. In this regard, the probe usage limitation component 51 is able to store an initial probe usage time T_i and a predetermined elapsed time limit T_L . The initial probe usage time T_i is written into the probe usage limitation component 51 by the generator 12 when the medical probe 302 is initially operated. Subsequently, the initial probe usage time T_i can be read from the probe usage limitation component 51 and used to limit re-usage of the medical probe 302, based on the initial operation of the medical probe 302.

Utility of storing the initial probe usage time T_i within the probe usage limitation component 51, however, is not to be limited to limiting re-usage of the medical probe 302, but may extend to any application wherein knowledge of the time of the initial operation of the medical probe 302 would be useful. For example, the initial probe usage time T_i can be read from a faulty medical probe to facilitate troubleshooting of the medical probe.

The predetermined elapsed time limit T_L is written into the probe usage limitation component 51, preferably during the manufacturing process. The predetermined elapsed time limit T_L represents the amount of time after the medical probe has been initially

operated that the medical probe remains operable. Preferably, the value of the predetermined elapsed time limit T_L is selected to minimize re-use of the medical probe 302, while still allowing legitimate use of the medical probe 302. For example, the value of the predetermined elapsed time limit T_L can be selected to be between the maximum
5 time expected to perform a therapeutic ablation procedure on a patient and the minimum time expected to complete a sterilization cycle on the medical probe 302. By way of non-limiting example, the value of the predetermined elapsed time limit T_L can be selected to be 24 hours.

In the illustrated embodiment, the probe usage limitation component 51 comprises
10 an EEPROM. Preferably, the probe usage limitation component 51 includes temperature-sensing capability to provide a known reference temperature T_r for the temperature sensor 19, which, in the illustrated embodiment, is a thermocouple having a first temperature-sensing element 21 located adjacent the ablation electrode 16, and a second temperature-sensing element 23 located in the handle 20 (shown in Fig. 1) adjacent the probe usage
15 limitation component 51. In this manner, the need for a thermistor, which was previously used to measure the reference temperature adjacent the second temperature-sensing element 23, is obviated. By way of non-limiting example, a Model DS2434 battery identification chip manufactured by Dallas Semiconductor is an EEPROM that provides temperature-sensing capability.

20 The structure of the probe usage limitation component 51 is not limited to an EEPROM, but can take the form of any non-volatile component that allows storage of data therein on-the-fly. Alternatively, the structure of the probe usage limitation component 51

can comprise other microchips, such as, e.g., non-volatile RAM or volatile RAM in conjunction with a battery.

The generator 12 includes an RF power source 52, which is connected through a transformer 54 to first and second conducting lines 56 and 58. In the illustrated

5 environment, the power source 52 delivers up to 150 watts of power at a frequency of 500kHz. The first conducting line 56 leads to the ablation electrode 16, and the second conducting line 58 leads to the indifferent electrode 18. The generator 12 further includes an interlock device 60, which is electrically coupled to the power source 52. Activation of the interlock device 60 via a control signal prevents the power source 52 from outputting
10 power to the transformer 54, thereby preventing delivery of RF energy to the tissue to be ablated.

The generator 12 includes a clock 62, which, as will be discussed in further detail below, generates a time signal T for time-keeping functions. In the illustrated embodiment, the time signal T represents an absolute time, such as, e.g., what one may obtain from a
15 standard household clock.

The generator 12 further includes a control circuit 68, which is electrically coupled to the clock 62 to obtain the time signal T. The control circuit 68 is electrically coupled to the probe usage limitation component 51 of the medical probe 302 when the medical probe 302 is physically connected to the generator 12, so that the control circuit 68 can
20 read and write data to and from the probe usage limitation component 51.

The control circuit 68 is configured for deriving an initial probe usage time T_i from the time signal T when the medical probe is initially operated, and for writing the initial probe usage time T_i to the probe usage limitation component 51, assuming that the probe

usage limitation component 51 has not previously been initialized with the initial probe usage time T_i . In this manner, the initial probe usage time T_i is reported to the generator 12 or a like generator if the medical probe 302 is reconnected to the generator 12, i.e., when the medical probe 302 is physically disconnected from the generator 12 and then physically reconnected to the generator 12 or like generator.

The control circuit 68 is further configured for reading the initial probe usage time T_i from the probe usage limitation component 51, and generating an interlock control signal S_c based on the initial probe usage time T_i . In particular, the control circuit 68 is configured for reading the predetermined elapsed time limit T_L from the probe usage limitation component 51 and determining if the predetermined elapsed time limit T_L has expired. The control circuit 68 determines expiration of the predetermined elapsed time limit T_L by determining an elapsed probe usage time T_E based on a comparison between initial probe usage time T_i and a reference time T_R obtained from the time signal T , i.e., the initial probe usage time T_i is subtracted from the reference time T_R . In the illustrated embodiment, the reference time T_R represents the current absolute time. The control circuit 68 then compares the elapsed probe usage time T_E with the predetermined elapsed time limit T_L to determine if the predetermined elapsed time limit T_L has expired, i.e., if the elapsed probe usage time T_E exceeds the predetermined elapsed time limit T_L , the predetermined elapsed time limit T_L is considered expired. For the purposes of this invention, the value of the elapsed probe usage time T_E exceeds the value of the predetermined elapsed time limit T_L if the value of the elapsed probe usage time T_E is equal to or greater than the value of the predetermined elapsed time limit T_L .

The control circuit 68 is electrically coupled to the interlock device 60 and applies the interlock control signal S_c thereto. Application of the interlock control signal S_c activates the interlock device 60, preventing power from being outputted from the power source 52, and in turn, preventing subsequent conveyance of RF energy to tissue in contact with the ablation electrode 16.

With reference to Fig. 3, operation of the medical probe system 300, in identifying and disabling re-used single use devices based on the time elapsed from the first therapeutic use of the medical probe 302, will now be described. Without knowledge of whether the medical probe 302 has been previously operated, the medical probe system 300 is set up by physically connecting the medical probe 302 to the generator 12 and powering on the generator 12 (step 70). At this point, the generator 12 is in standby mode.

The generator 12 then reads the data from the probe usage limitation component 51 (step 72). After insertion of the medical probe 302 within the patient's body and placement of the ablation electrode 16 in contact with the tissue to be ablated, the RF power control button 44 is depressed to conditionally initiate power transmission from the generator 12 and subsequent delivery of RF energy to the tissue (step 74). At this point, the generator 12 is in delivery mode.

The generator 12 then determines whether the probe usage limitation component 51 has already been initiated with an initial probe usage time T_i by analyzing the data read from the probe usage limitation component 51 (step 76). If the probe usage limitation component 51 has not already been initialized with the initial probe usage time T_i , the generator 12 initializes the probe usage limitation component 51 with the initial probe usage time T_i (step 78). That is, the generator 12 derives the initial probe usage time T_i .

from the timing signal T , and writes the initial probe usage time T_i into the probe usage limitation component 51 of the medical probe 302. In this manner, if the medical probe 302 is physically disconnected and reconnected to the generator 12 or like generator and conditionally operated, the generator 12 or like generator can ascertain when the medical probe 302 was initially operated. Since the medical probe 302 has not been previously operated, the generator 12 then allows operation of the medical probe 302, i.e., the interlock control signal S_c is not generated (step 88).

In the illustrated embodiment, initialization of the probe usage limitation component 51 occurs upon effective operation of the medical probe 302, i.e., operation of the medical probe 302 in such a manner as to form a lesion on the tissue. In this regard, to prevent premature initiation of the probe usage limitation component 51, the generator 12 does not write to the probe usage limitation component 51 during non-therapeutic operation of the medical probe 302, i.e., when operation of the medical probe 302 does not result in a tissue lesion. Such probe operation can be caused by a variety of reasons, including inadequate contact between the ablation electrode 16 and the tissue, inadequate energy delivery to the ablation electrode 16 and inadequate duration of energy application.

After initiation of the probe usage limitation component 51 with the initial probe usage time T_i (step 78), or if the generator 12 determines that the probe usage limitation component 51 has already been initialized with the initial probe usage time T_i (step 76), the generator 12 determines whether the predetermined elapsed time limit T_L has expired (steps 80-84). First, the reference time T_R is determined from the time signal T (step 80). Then, the elapsed probe usage time T_E is determined by subtracting the initial probe usage

time T_i from the reference time T_R (step 82). Then, the elapsed probe usage time T_E is compared to the predetermined elapsed time limit T_L (step 84).

If the elapsed probe usage time T_E exceeds the predetermined elapsed time limit T_L , the generator 12 prevents operation of the medical probe 302, i.e., the interlock control signal S_C is generated, which is subsequently transmitted to the interlock device 60 (step 5 86). The medical probe 302 is thus rendered inoperable and cannot be re-used, thereby preventing further operation of the medical probe 302. If the value of the elapsed probe usage time T_E does not exceed the value of the predetermined elapsed time limit T_L , the generator 12 allows operation of the medical probe 302, i.e., the interlock control signal S_C 10 is not generated (step 88).

The generator 12 periodically (e.g., every second) determines if the elapsed probe usage time T_E has expired and conditionally allows operation of the medical probe 302, based on the expiration of the elapsed probe usage time T_E , until the medical probe 302 is disconnected from the generator 12 (steps 80-88). Alternatively, the expiration of the 15 elapsed probe usage time T_E is determined only one time per probe connection. That is, once the elapsed probe usage time T_E is determined not to be expired, the medical probe 302 can be operated without limitation until the medical probe 302 is physically disconnected from, and again connected, to the generator 12 (step 90).

In this manner, the physician has a certain amount of time in which to operate the 20 medical probe 302 after initiation of the probe usage limitation component 51, notwithstanding the fact that the medical probe 302 can be repeatedly disconnected from and reconnected to the generator 12. Thus, re-usage of the medical probe 302 can be minimized without prematurely rendering the medical probe 302 inoperable.

Limiting Re-Use of Devices Based On Detecting Environment Changes

Referring to Fig. 4, the medical probe system 300 is capable of identifying and disabling re-used single use devices based on detecting environmental changes. In this regard, the probe usage limitation component 51 is capable of detecting an environmental condition to which the medical probe 302 is exposed, and determines if the environmental condition is indicative of exposure of the medical probe 302 to a sterilization cycle.

In a typical sterilization procedure, the medical probe 302 is placed inside a breathable pouch, and a biological indicator is applied to the medical probe 302 according to a validated hospital procedure. The medical probe 302 and pouch are then placed inside a sterilization chamber and exposed to a sterilization cycle. The sterilization cycle may, for example, include preconditioning, exposure, post-vacuum, and aeration steps, which are performed under a variety of environmental conditions. For example, the preconditioning step may involve subjecting the medical probe 302 to a temperature of 125-145°F, humidity of 55-75%, and a pressure (vacuum) at 1.9-3.9 p.s.i.a. for 30-45 minutes. The exposure step may involve subjecting the medical probe 302 to an oxide, such as, e.g., 100% ethylene oxide (EO) at 600 mg/L, and a temperature of 125-145°F for four hours. The post-vacuum step may involve twice subjecting the medical probe 302 to a pressure (vacuum) of 1.9-3.9 p.s.i.a. The aeration step may involve subjecting the medical probe 302 to a temperature of 120°F-145°F for twelve hours.

Thus, if the probe usage limitation component 51 detects an environmental condition that is indicative of a sterilization cycle, such as, e.g., temperature, humidity, pressure or chemical, then the probe usage limitation component 51 stores a probe

sterilization indicator I_s indicating that the medical probe 302 has been sterilized. For example, if the probe usage limitation component 51 detects one or more of the following conditions: a temperature above 125°F, a humidity above 55%, a pressure below 3.9 p.s.i.a., or an oxide, such as, e.g., EO, the probe usage limitation component 51 stores a probe sterilization indicator I_s .

Preferably, as part of the post-manufacturing process, the medical probe 302 is subjected to an initial sterilization cycle. In this case, the medical probe 302 will be shipped to the physician with the probe sterilization indicator I_s stored in the probe usage limitation component 51. The control circuit ultimately limits re-usage of the medical probe 302, based on a detection of a re-sterilization of the medical probe 302, an indication that the physician is attempting to operate the medical probe 302 illegitimately.

A medical probe 302 that has been re-sterilized (the operation of which would be considered illegitimate operation of the medical probe 302) must be distinguished from a medical probe 302 that has been sterilized only once as a part of the post-manufacture process (the operation of which would be considered legitimate). In this regard, as will be discussed in further detail below, the probe usage limitation component 51 is configured for storing a previous probe usage indicator I_U , indicating that the medical probe 302 has been previously operated and allowing post-manufacturing sterilization to be distinguished from re-sterilization of the medical probe 302.

In the illustrated embodiment, the probe usage limitation component 51 comprises battery-operated non-volatile RAM with environmental sensing capabilities. A non-limiting example of such a probe usage limitation component 51 is an electrolytic sensor or chemFET.

The control circuit 68 is configured for reading data from the probe usage limitation component 51 and determining if the data comprises the probe sterilization indicator I_s and the previous probe usage indicator I_U , and generating an interlock control signal S_C based on the presence of the probe sterilization indicator I_s and the previous probe usage indicator I_U within the data read from the probe usage limitation component 51. In particular, the interlock control signal S_C is generated if both the probe sterilization indicator I_s and the previous probe usage indicator I_U are present, or if both the probe sterilization indicator I_s and the previous probe usage indicator I_U are not present. The interlock control signal S_C is not generated if the probe sterilization indicator I_s is present and the previous probe usage indicator I_U is not present, or if the probe sterilization indicator I_s is not present and the previous probe usage indicator I_U is present.

With reference to Fig. 5, operation of the medical probe 302 during sterilization will now be described. After manufacture of the medical probe 302, the probe usage limitation component 51 carried by the medical probe 302 is continuously detecting an environmental condition to which the medical probe 302 is exposed (step 92). The environmental condition detected by the probe usage limitation component 51 can be, e.g., temperature, pressure, moisture and/or chemical-based. The probe usage limitation component 51 then determines if the environmental condition detected is associated with a sterilization cycle, i.e., if the environmental condition matches a known environmental condition to which probes are subjected during the sterilization cycle (step 94). If the probe usage limitation component 51 detects an environmental condition associated with the sterilization cycle, such as, e.g., a temperature of 130°F, the probe usage limitation component 51 stores the probe sterilization indicator I_s , indicating that the medical probe

302 has been either sterilized, whether it be legitimately during the post-manufacture process or illegitimately during a re-sterilization process (step 96). This sterilization indicator I_s is reported later to the generator 12 to indicate that the medical probe 302 has been sterilized. The probe usage limitation component 51 continuously detects and
5 determines if the environmental condition to which the medical probe 302 is exposed is associated with a sterilization cycle.

With reference to Fig. 6, operation of the medical probe system 300, in identifying and disabling re-used single use devices based on detecting environmental changes, will now be described. Without knowledge of whether the medical probe 302 has been
10 sterilized or previously operated, the medical probe system 300 is set up by physically connecting the medical probe 302 to the generator 12 and powering on the generator 12 (step 98). At this point, the generator 12 is in standby mode. The generator 12 then reads the data from the probe usage limitation component 51 (step 100). After insertion of the medical probe 302 within the patient's body and placement of the ablation electrode 16 in
15 contact with the tissue to be ablated, the RF power control button 44 (shown in Fig. 1) is depressed to conditionally initiate power transmission from the generator 12 and subsequent delivery of RF energy to the tissue (step 102). At this point, the generator 12 is in delivery mode.

The generator 12 then determines whether the probe sterilization indicator I_s is
20 present in or absent from the data read from the probe usage limitation component 51, and thus whether the medical probe 302 has been sterilized (step 104). If the probe sterilization indicator I_s is present, the generator 12 determines whether the previous probe usage indicator I_u is present in the data read from the probe usage limitation component

51, and thus whether the medical probe 302 has been previously operated (step 106). If the previous probe usage indicator I_U is absent, the medical probe 302 is determined to be sterilized and unused. Thus, the generator 12 allows operation of the medical probe 302, i.e., the interlock control signal S_C is not generated (step 108).

5 The generator 12 then writes the previous probe usage indicator I_U to the probe usage limitation component 51 (step 110). In this manner, if the medical probe 302 is physically disconnected and reconnected to the generator 12 or like generator and conditionally operated, either legitimately or illegitimately, the generator 12 or like generator can ascertain that the medical probe 302 has been previously operated. The
10 generator 12 also clears the probe sterilization indicator I_S from the probe usage limitation component 51 (step 112). In this manner, if the medical probe 302 is physically disconnected and reconnected to the generator 12 without undergoing a sterilization cycle, the medical probe 302 will not be prematurely rendered inoperable, which would otherwise result from the presence of the previous probe usage indicator I_U . That is, as will
15 be discussed in further detail below, if the probe sterilization indicator I_S is absent and the previous probe usage indicator I_U is present, the medical probe 302 can still be operated.

In the illustrated embodiment, storage of the previous probe usage indicator I_U into the probe usage limitation component 51 occurs upon effective operation of the medical probe 302, i.e., operation of the medical probe 302 in such a manner as to form a lesion
20 on the tissue. In this regard, to prevent premature storage of the previous probe usage indicator I_U into the probe usage limitation component 51, the generator 12 does not write to the probe usage limitation component 51 during non-therapeutic operation of the medical probe 302, i.e., when faulty operation of the medical probe 302 does not result in

a tissue lesion. Such faulty probe operation can be caused by a variety of reasons including, inadequate contact between the ablation electrode 16 and the tissue, inadequate energy delivery to the ablation electrode 16 and inadequate duration of energy application.

If the previous probe usage indicator I_U is present (step 106), the medical probe 302 is determined to have been previously used and re-sterilized illegitimately, and thus, the generator 12 prevents operation of the medical probe 302, i.e., the interlock control signal S_C is generated (step 114).

If the probe sterilization indicator I_S is absent (step 104), the generator 12 determines if the previous probe usage indicator I_U is present in, or absent from, the data read from the probe usage limitation component 51 (step 116). If the previous probe usage indicator I_U is present, the medical probe 302 is determined to have been physically disconnected from, and reconnected to, the generator 12 in the middle of a procedure, and thus further operation of the medical probe 302 is legitimate. As such, the generator 12 allows operation of the medical probe 302; i.e., the interlock control signal S_C is not generated (step 118).

If, after the probe sterilization indicator I_S is determined to be absent, the previous probe usage indicator I_U is absent, the medical probe 302 is determined to be improperly initially sterilized or not sterilized at all, and thus, the generator 12 prevents operation of the medical probe 302, i.e., the interlock control signal S_C is generated (step 120).

In this manner, as long as the medical probe 302 is not re-sterilized, the physician can operate the medical probe 302 after initiation of the probe usage limitation component 51, notwithstanding the fact that the medical probe 302 can be repeatedly disconnected

from and reconnected to the generator 12. Thus, re-usage of the medical probe 302 can be minimized without prematurely rendering the medical probe 302 inoperable.

The medical probe system 300 can use other probe indicators in combination with the probe sterilization indicator I_s to identify and disable re-used single use devices based on detecting environmental changes. For example, when the medical probe 302 is manufactured, the estimated sterilization date of the medical probe 302, or alternatively, the manufacture date of the medical probe 302, can be stored in the probe usage limitation component 51. A typical estimated probe sterilization date is a month after the manufacture date of the probe. Upon actual sterilization of the medical probe 302, the probe sterilization indicator I_s and the date of actual sterilization are stored in the probe usage limitation component 51. Upon connection of the medical probe 302 to the RF generator 12, if the probe sterilization indicator I_s is absent, the medical probe 302 is determined to be un-sterilized, and therefore is prevented from being operated. If the sterilization indicator I_s is present and the actual date of sterilization is later the estimated actual date of sterilization, the medical probe 302 is determined to be re-sterilized, and therefore is prevented from being operating.

The probe sterilization indicator I_s cannot only be used to limit re-usage of medical probes 302, but can also be used to enforce the shelf-life of medical probes 302. In this regard, upon sterilization of the medical probe 302, the sterilization indicator I_s and actual date of sterilization can be stored in the probe usage limitation component 51. If, upon connection of the medical probe 302 to the RF generator 12, the difference between the actual sterilization date and the current date exceeds a predetermined period of time (i.e.,

the shelf-life of the medical probe 302), the medical probe 302 is prevented from being operated.

Limiting Re-Use of Devices Based On Cataloging Catheter Usage

5 Referring to Figs. 7 and 8, the medical probe system 300 is alternatively capable of identifying and disabling re-used single use devices based on cataloging catheter usage. In this regard, the control circuit 68 includes a catalog 122 (shown in further detail in Fig. 8), obviating the need for a probe usage limitation component 51 within the medical probe. In this embodiment, the control circuit 68 is configured for storing historical operational
10 indicators I_O obtained from a multitude of medical probes 302 previously operated with the generator 12. In the illustrated embodiment, a set of historical operational indicators I_O are stored for each of the multitude of previously operated medical probes 302. Each set of historical operational indicators I_O includes a probe identification code C_i and supplementary historical operational indicators I_{SO} . Each probe identification code C_i not
15 only provides historical operational information about the respective medical probe 302, but also provides a means to cross-reference the supplementary historical operational indicators I_{SO} within the respective set of historical operational indicators I_O , as will be described in further detail below.

The supplementary operational indicators I_{SO} provide historical operational
20 information concerning the operation of a medical probe 302, in addition to that provided by the presence of a cataloged probe identification code C_i . By way of non-limiting example, a supplementary operational indicator I_O can be the initial probe usage time T_i described previously above. Alternatively, a supplementary operational indicator I_{SO} can

be an incremental probe usage U_i , which represents the number of times a respective medical probe 302 has been operated. As will be described in further detail below, the historical operational indicators I_{SO} can provide historical operational information of a selected medical probe 302 to a physician. The historical operational indicators I_{SO} can be further used to conditionally operate the selected medical probe 302.

The catalog 122 is configured for retaining the probe historical operational indicators I_O , notwithstanding disconnection of the medical probe 302 from the generator 12. Thus, storage of the historical operational indicators I_O within the catalog 122 for subsequent use obviates the need to store historical operational information within the medical probe 302 itself, thereby minimizing probe expense. For instance, absent the requirement to store historical operational information in the electronic storage componentry of the medical probe 302, a ROM or one or more resistors, which are relatively inexpensive, can be used as the probe identification component 50, merely for storage of the probe identification code C_i .

The control circuit 68 is configured for, upon connection of a selected medical probe 302 to the generator 12, reading the probe identification code C_i from the probe identification component 50 of the selected medical probe 302. The control circuit 68 is further configured for, upon initial operation of the selected medical probe 302, storing the probe identification code C_i in the catalog 122 as a historical operational indicator I_O , the presence of which indicates that the selected medical probe 302 has been previously operated with the generator 12.

The control circuit 68 is also configured for obtaining additional historical operational information concerning the operation of the selected medical probe 302, and

storing this information in the catalog 122 as supplemental historical operational indicators

15 I_{so} . For instance, upon initial operation of the selected medical probe 302, the control circuit 68 can obtain the time of initial operation of the selected medical probe 302, which can then be stored in the catalog 122 as an initial probe usage time T_i next to the probe identification code C_i corresponding to the selected medical probe 302. Also, each time the selected medical probe 302 is operated, the control circuit 68 can obtain a cumulative number of times the selected medical probe 302 has been operated, which can then be stored in the catalog 122 as an incremental probe usage U_i next to the probe identification code C_i corresponding to the selected medical probe 302. Thus, a set of historical operational indicators I_o for the selected medical probe 302 can be stored in the catalog 122.

The control circuit 68 is also configured for recalling historical operational information about a selected medical probe 302, which has been stored in the catalog 122 as the set of historical operational indicators I_o . In particular, the control circuit 68 can determine if the medical probe 302 has been previously operated with the generator 12 by determining if the probe identification code C_i read from the medical probe 302 matches any of the probe identification codes C_i stored in the catalog 122. The control circuit 68 can also determine additional historical operational information, such as, e.g., the initial time of operation of the medical probe 302 or the number of times the medical probe 302 has been previously operated, by obtaining the supplemental operational indicators I_o , associated with the matched probe identification code C_i . The recalled historical operational information can optionally be displayed on the display (not shown) for viewing by the physician.

The control circuit 68 is further configured for conditionally operating the selected medical probe 302, based on the set of historical operational indicators I_O corresponding to the selected medical probe 302. That is, the control circuit 68 is configured for generating an interlock control signal S_C based on the probe identification code C_i and/or associated supplemental operational indicators I_{SO} . As will be discussed in further detail below, conditional operation of the medical probe 302 can not only be accomplished during the initial connection between the medical probe 302 and the generator 12, but can also be accomplished during reconnection between the medical probe 302 and the generator 12 (by recalling from the catalog 122), since the cataloged historical operational information is not lost upon disconnection of the medical probe 302 from the generator 12.

By way of non-limiting example, the control circuit 68 can conditionally operate the selected medical probe 302, based on the probe identification code C_i and the probe sterilization indicator I_s , which is stored in the selected medical probe 302 and read out by the control circuit 68, as previously discussed above. In particular, as discussed above, a presence of the probe identification code C_i in the catalog 122 indicates that the selected medical probe 302 has been previously operated. A presence of the probe sterilization indicator I_s in the selected medical probe 302 indicates that the selected medical probe 302 has been previously sterilized. In this manner, use of the probe identification code C_i allows a medical probe that has been legitimately sterilized by a post-manufacture process to be distinguished from a medical probe that has been illegitimately re-sterilized after usage. Upon operation of the selected medical probe 302, the control circuit 68 stores the probe identification code C_i in the catalog 122 and clears the probe sterilization indicator I_s from the selected medical probe 302, so that the selected medical probe 302 is not

prematurely rendered inoperable through legitimate adjustment of the connection between the selected medical probe 302 and the generator 12.

Operation of the selected medical probe 302 is prevented if both the probe sterilization indicator I_s and the probe identification code C_i are present (indicating that the selected medical probe 302 has been re-sterilized), or if both the probe sterilization indicator I_s and the probe identification code C_i are not present (indicating that the selected medical probe 302 has never been sterilized). Operation of the selected medical probe 302 is allowed if the probe identification code C_i is present and the probe identification code C_i is not present (indicating that the medical probe has only been sterilized once and never operated), or the probe sterilization indicator I_s is not present and the probe identification code C_i is present (indicating that the medical probe has only been sterilized once and is currently being used in an initial procedure).

By way of further non-limiting example, the control circuit 68 can conditionally operate the selected medical probe 302, based on the initial probe usage time T_i and the predetermined elapsed time limit T_L , which may be encoded in the probe identification code C_i of the selected medical probe 302, as previously discussed above. Upon initial operation of the selected medical probe 302, the control circuit 68 stores the initial probe usage time T_i , along with the probe identification code C_i , in the catalog 122. The control circuit 68 prevents operation of the selected medical probe 302 when the predetermined elapsed time limit T_L is expired. Expiration of the predetermined elapsed time limit T_L occurs when the difference between the initial probe usage time T_i and a reference time exceeds the predetermined elapsed time limit T_L .

By way of further non-limiting example, the control circuit 68 can conditionally operate the selected medical probe 302, based on the incremental probe usage U_i and a predetermined maximum usage limit U_L , which may be encoded in the probe identification code C_i of the selected medical probe 302. Upon initial operation of the selected medical probe 302, the control circuit 68 stores a "1" as the incremental probe usage U_i , along with the probe identification code C_i , in the catalog 122. Upon each subsequent operation of the selected medical probe 302, the control circuit 68 increments the incremental probe usage U_i by one. The control circuit 68 prevents operation of the selected medical probe 302 when the incremental probe usage U_i exceeds the predetermined maximum usage limit U_L .

Further details regarding the use of an incremental probe usage U_i and a predetermined maximum usage limit U_L is further disclosed in U.S. Patent No. 5,383,874, entitled "Systems for Identifying Catheters and Monitoring Their Use," which issued January 24, 1995, and which is hereby fully and expressly incorporated herein by reference.

With reference to Fig. 9, operation of the medical probe system 300, in identifying and disabling re-used single use devices based on cataloging catheter usage, will now be described. The medical probe system 300 is set up by physically connecting a selected medical probe 302 (the historical operation of which is not known) to the generator 12 and powering on the generator 12 (step 124). At this point, the generator 12 is in standby mode. The generator 12 then reads the probe identification code C_i from the probe identification component 50 of the selected medical probe 302, i.e., obtains the selected probe identification code C_i (step 126). After insertion of the selected medical probe 302

within the patient's body and placement of the ablation electrode 16 in contact with the tissue to be ablated, the RF power control button 44 (shown in Fig. 1) is depressed to conditionally initiate power transmission from the generator 12 and subsequent delivery of RF energy to the tissue (step 128). At this point, the generator 12 is in delivery mode.

5 The generator 12 then obtains historical operational information for the selected medical probe 302 from the catalog 122 (steps 130-134). That is, the catalog 122 is searched to determine if the selected probe identification code C_i matches a cataloged probe identification code C_i (step 130). If the selected probe identification code C_i does match a cataloged probe identification code C_i , the set of historical operational indicators I_o associated with the matched identification code C_i (i.e., the set of historical operational indicators I_o corresponding to the selected medical probe 302) is obtained (step 132).

As discussed above, the set of historical operational indicators I_o comprises the probe identification code C_i , the presence of which indicates that the selected medical probe 302 has been previously operated; the initial probe usage time T_i indicating that time of initial operation of the selected medical probe 302; and the incremental probe usage U_i , indicating the number of times the selected medical probe 302 has been previously operated. This historical operational information can be optionally displayed on the display (not shown) for viewing by the physician.

15 If the selected probe identification code C_i does not match a cataloged probe identification code C_i , the selected medical probe 302 is deemed to have not been previously operated, and a set of historical operational indicators I_o corresponding to the selected medical probe 302 does not exist. This historical operational information can be optionally displayed on the display (not shown) for viewing by the physician.

The generator 12 then determines whether or not the selected medical probe 302 should be operated, based on the historical operational information obtained from the catalog 122 (step 134). For example, as discussed above, the generator 12 can conditionally operate the selected medical probe 302, based on the presence of the probe identification code C_i in the catalog 122 and a presence of the probe sterilization indicator I_s in the selected medical probe 302; and/or based on the initial probe usage time T_i and the predetermined elapsed usage limit T_L ; and/or the incremental probe usage U_i and the predetermined maximum usage limit U_L .

If the generator 12 determines that the selected medical probe 302 can be operated, the generator 12 allows operation of the selected medical probe 302, i.e., the interlock control signal S_c is not generated (step 136). If operation of the selected medical probe 302 is the first time that the selected medical probe 302 has been operated, i.e., if the selected probe identification code C_i did not match a cataloged probe identification code C_l (step 138), the catalog 122 is initialized with a set of historical operational indicators I_o corresponding to the selected medical probe 302, based on the operation of the selected medical probe (step 140). That is, the selected probe identification code C_i , along with the time of initial operation and number of times the selected medical probe 302 has been operated (1 time in this case), is stored in the catalog 122.

In the illustrated embodiment, initialization of the catalog 122 with the set of historical operational indicators I_o occurs upon effective operation of the medical probe 302, i.e., operation of the medical probe 302 in such a manner as to form a lesion on the tissue. In this regard, to prevent premature initialization of the catalog 122, the set of historical operational indicators I_o is not stored in the catalog 122 during non-therapeutic

operation of the medical probe 302, i.e., when faulty operation of the medical probe 302 does not result in a tissue lesion. Such faulty probe operation can be caused by a variety of reasons, including inadequate contact between the ablation electrode 16 and the tissue, inadequate energy delivery to the ablation electrode 16 and inadequate duration of energy application.

If operation of the selected medical probe 302 is not the first time that the selected medical probe 302 has been operated, i.e., if the selected probe identification code C_i did match a cataloged probe identification code C_i (step 138), the generator 12 updates the set of historical operational indicators I_o corresponding to the selected medical probe 302 (step 142). For instance, the incremental probe usage U_i may be incremented by one. If the generator 12 determines that the selected medical probe 302 cannot be operated (step 134), the generator 12 prevents operation of the selected medical probe 302, i.e., the interlock control signal S_c is generated (step 144).

The generator 12 periodically determines (e.g., every second, when operation of the selected medical probe 302 is based on the initial probe usage time T_i , and every time the selected medical probe 302, when operation of the selected medical probe 302 is based on the incremental probe usage U_i) whether the selected medical probe 302 should be operated, and conditionally allows operation of the medical probe 302, based on the set of historical operational indicators I_o , until the selected medical probe 302 is disconnected from the generator 12 (steps 134-144). Alternatively, the generator 12 determines whether the selected medical probe 302 should be operated only one time per probe connection. That is, once operation of the selected medical probe 302 is allowed, the selected medical

probe 302 can be operated without limitation until the selected medical probe 302 is physically disconnected from, and again connected to, the generator 12 (step 146).

Preventing Automatic Identification of Devices

5 Referring to Fig. 10, the medical probe system 300 is alternatively capable of preventing automatic identification of re-used single use devices. In this regard, the probe usage limitation component 51 is a disabling circuit, the activation of which disables the probe identification component 50 by modifying the manner in which the probe identification component 50 functions. In the illustrated embodiment, the disabling circuit 10 51 is a switch that changes the output of the probe identification component 50. For example, if the probe identification component 50 comprises one or more resistors, the disabling circuit 51 can change the output resistance value seen by the generator 12 by shorting the output of the one or more resistors. Of course, if probe identification component 50 comprises a microchip, one or more pins of the microchip can be shorted 15 to create the same effect.

Disabling of the probe identification component 50 prevents the correct probe identification code C_i from being read therefrom, which, as will be discussed in further detail below, precludes the medical probe 302 from being operated. Prior to initial operation of the medical probe 302, the disabling circuit 51 is configured to enable the 20 probe identification component 50. That is, the disabling circuit 51 provides an open circuit to the output of the probe identification component 50. A non-limiting example of a disabling circuit 51 is an electrolytic switch. Application of a disable signal S_d on the disabling circuit 51 from an external source, i.e., the generator 12, activates the disabling

circuit 51. Preferably, once the disabling circuit 51 is activated, disablement of the probe identification component 50 is continued even when the power is removed from the disabling circuit 51. In this manner, no generator 12 will be able to read the probe identification code C_i from the probe identification component 50 after the probe
5 identification component 50 has been disabled.

The control circuit 68 is configured for determining whether the medical probe 302 has been previously used, i.e., whether the medical probe 302 has been used to render a complete procedure on a patient. The control circuit 68 can determine whether the medical probe 302 has been previously used by employing a variety of techniques.

10 By way of non-limiting example, the control circuit 68 can base the determination of whether the medical probe 302 has been previously used on the time elapsed after an initial operation of the medical probe 302, as discussed above. In particular, upon initial operation of the medical probe 302, the control circuit 68 stores an initial probe usage time T_i in the medical probe 302 (either in the probe identification component 50 or other
15 electronic storage componentry), or in a catalog within the control circuit 68, along with the probe identification code C_i , which can then be later read out by the control circuit 68.

The control circuit 68 then determines if a predetermined elapsed time limit T_L (which may be encoded into the probe identification code C_i) has expired, i.e., if the difference between the initial probe usage time T_i and a reference time exceeds the predetermined
20 elapsed time limit T_L . The medical probe 302 is determined to have been previously used when the predetermined elapsed time limit T_L has expired.

By way of further non-limiting example, the control circuit 68 can base the determination of whether the medical probe 302 has been previously used on a number of

times the medical probe 302 has been operated, as discussed previously. In particular, upon initial operation of the medical probe 302, an incremental probe usage U_i can be stored in the medical probe 302 (either in the probe identification component 50 or other electronic storage componentry), or a catalog, along with the probe identification code C_i , which can then be later read out by the control circuit 68. Each time the medical probe 302 is operated, i.e., each time power is applied to the medical probe 302, the control circuit 68 increments the incremental probe usage U_i by one. The control circuit 68 then determines if the incremental probe usage U_i exceeds a predetermined probe usage limit U_L (which may be encoded into the probe identification code C_i). The medical probe 302 is determined to have been previously used when incremental probe usage U_i exceeds the predetermined probe usage limit U_L .

By way of further non-limiting example, the control circuit 68 can base the determination of whether the medical probe 302 has been previously used on both a previous operation of the medical probe 302 and a previous sterilization of the medical probe 302, as discussed previously. In particular, upon initial operation of the medical probe 302, a probe usage indicator I_U is stored in the medical probe 302 (either in the probe identification component 50 or other electronic storage componentry), or a catalog in the form of the probe identification code C_i , which can then be later read out by the control circuit 68. Upon sterilization of the medical probe 302, a probe sterilization indicator I_S is stored in the medical probe 302, preferably in an environmentally sensing microchip, which can be later read out by the control circuit 68. In this manner, use of the probe identification code C_i allows a medical probe that has been legitimately sterilized by

a post-manufacture process to be distinguished from a medical probe that has been illegitimately re-sterilized after usage.

Upon operation of the medical probe 302, the control circuit 68 stores the probe identification code C_i in the medical probe 302 or catalog 122 and clears the probe sterilization indicator I_s from the medical probe 302 or catalog, so that the medical probe 302 is not prematurely rendered inoperable through legitimate adjustment of the connection between the medical probe 302 and the generator 12. The medical probe 302 is determined to have been previously used if both the probe sterilization indicator I_s and the previous probe usage indicator I_u are present (indicating that the medical probe 302 has been re-sterilized), or if both the probe sterilization indicator I_s and the previous probe usage indicator I_u are not present (indicating that the medical probe 302 has never been sterilized).

As will be discussed in further detail below, previous probe use can be periodically determined after connection of the medical probe 302 to the generator 12, but before disconnection of the medical probe 302 from the generator 12. Alternatively, previous probe use is only determined once upon connection of the medical probe 302 to the generator 12.

The control circuit 68 is further configured for generating the disable signal S_D based on the determination of whether the medical probe 302 has been previously used. In particular, if the medical probe 302 is determined to have been previously used, the control circuit 68 generates the disable signal S_D , thereby activating the disabling circuit 51. As discussed above, activation of the disabling circuit 51 disables the probe identification component 50, thereby preventing the probe identification code C_i from

being read from the probe identification component 50. If the medical probe 302 is determined not to have been previously used, the control circuit 68 does not generate the disable signal S_D , thereby maintaining the inactivation of the disabling circuit 51. As discussed above, inactivation of the disabling circuit 51 enables the probe identification component 50, thereby allowing the correct probe identification code C_i to be read from the probe identification component 50.

In alternative embodiments, the disabling circuit 51 can be eliminated by configuring the control circuit 68 to directly disable the probe identification component 50 by modifying the probe identification code C_i in the probe identification component 50 (shown by dashed line). Preferably, this can be accomplished by clearing the probe identification code C_i from the probe identification component 50. In this case, the probe identification component 50 preferably comprises a microchip, such as, e.g., an EEPROM, which is non-volatile, yet easily allows data to be written thereto. Alternatively, the probe identification component 50 can be a digital circuit that outputs a pulse train or a laser-trimmed clock. Even more alternatively, the probe identification component 50 can be a solid-state component, such as, e.g., a variable resistor, an array of resistors or a laser-trimmed resistor.

In further alternative embodiments, the disabling circuit 51 comprises a microprocessor that internally determines whether the medical probe 302 has been previously used (i.e., exclusive of a disable signal S_D transmitted from the generator 12), and disables the probe identification component 50, based on this determination. The probe-based microprocessor can disable the probe identification component 50 by modifying the functioning of the probe identification component 50 itself, such as by

shorting the probe identification component 50, or by modifying the probe identification code C_i storage by the probe identification component 50, such as by clearing the probe identification code C_i therefrom. In this manner, no additional processing need take place in the generator 12. Thus, a standard generator having the capability of reading a probe identification code from a medical probe and conditionally operating the medical probe based on the value of the probe identification code, can be used with the medical probe 302.

The probe-based microprocessor can sense power flow from the generator 12 during operation of the medical probe 302, thereby allowing the microprocessor to detect such probe usage parameters as previous probe operation, an initial time of probe operation, and a number of times of probe operation. The probe-based microprocessor may also have the capability of detecting sterilization of the medical probe 302. Thus, in this manner, the probe-based microprocessor can base the determination of whether the medical probe 302 has been previously operated on an elapsed time after an initial operation of the medical probe 302, on a number of times the medical probe 302 has been operated, or on a previous operation of the medical probe 302 and previous sterilization of the medical probe 302.

The control circuit 68 is further configured for reading the probe identification code C_i from the probe identification component 50 when enabled. The control circuit 68 compares the probe identification code C_i read from the probe identification component 50 with a look-up table of approved probe identification codes. If the probe identification code C_i does not match any of the approved probe identification codes in the look-up table (either because the wrong medical probe is being used with the generator 12, or the probe

identification code C_i read by the control circuit 68 has been altered due to previous probe use), the control circuit 68 does not allow operation of the medical probe 302. That is, the control circuit 68 generates an interlock control signal S_c .

If the probe identification code C_i does match any of the approved probe
5 identification codes in the look-up table (because the right medical probe is used and the probe identification code C_i read by the control circuit 68 has not been altered due to no previous probe use), the control circuit 68 allows operation of the medical probe 302. That is, the control circuit 68 does not generate an interlock control signal S_c . Further details regarding the conditional operation of a medical probe based on an automatic
10 reading of a probe identification code is further disclosed in U.S. Patent No. 5,383,874, which has previously been incorporated herein by reference.

The control circuit 68 is electrically coupled to the interlock device 60 and applies the interlock control signal S_c thereto. Application of the interlock control signal S_c activates the interlock device 60, preventing power from being outputted from the power
15 source 52, and in turn, preventing subsequent conveyance of RF energy to tissue in contact with the ablation electrode 16.

With reference to Fig. 11, operation of the medical probe system 300 will now be described. The medical probe system 300 is set up by physically connecting a medical probe 302 (the historical operation of which is not known) to the generator 12 and
20 powering on the generator 12 (step 148). At this point, the generator 12 is in standby mode. The generator 12 then reads the probe identification code C_i from the probe identification component 50 of the medical probe 302 (step 150). After insertion of the medical probe 302 within the patient's body and placement of the ablation electrode 16 in

contact with the tissue to be ablated, the RF power control button 44 (shown in Fig. 1) is depressed to conditionally initiate power transmission from the generator 12 and subsequent delivery of RF energy to the tissue (step 152). At this point, the generator 12 is in delivery mode.

5 The generator 12 then determines if the probe identification code C_i read from the probe identification component 50 matches an approved probe identification code located in the look-up list (step 154). If the probe identification code C_i does not match an approved probe identification code (either because the medical probe 302 is not approved for use with the generator 12 or the probe identification component 50 located in the
10 medical probe 302 has been disabled due to detection of previous probe usage), the generator 12 prevents operation of the medical probe 302, i.e., the interlock control signal S_c is generated (step 156). It should be noted that, if the probe identification component 50 is disabled by shorting, or if the probe identification code C_i is cleared from the probe identification component 50, the generator 12 will not be able to obtain a value from the
15 probe identification component 50, in which case the generator 12 need not access the look-up list to ultimately determine that the medical probe 302 should not be operated.

 If the probe identification code C_i does match an approved probe identification code (because the medical probe 302 is approved for use with the generator 12 and had not been determined to have been previously used during a previous connection), the
20 generator 12 determines if the medical probe 302 has been previously used based on, e.g., the presence of the probe identification code C_i and the presence of the probe sterilization indicator I_s ; and/or the initial probe usage time T_i and the predetermined elapsed usage

limit T_L ; or the incremental probe usage U_i and the predetermined maximum usage limit U_L , as discussed above (step 158).

If the medical probe 302 is determined not to have been previously used, the generator 12 allows operation of the medical probe 302, i.e., the interlock control signal S_C is not generated (step 160). If the medical probe 302 is determined to have been previously used, the generator 12 prevents operation of the medical probe 302, i.e., the interlock control signal S_C is generated (step 162). In addition, the generator 12 disables the probe identification component 50 in the medical probe 302 by transmitting the disable signal S_D to the disabling circuit 51, or alternatively, by clearing the probe identification code C_i from the probe identification component 50 (step 164). In this manner, when the medical probe 302 is reconnected to any generator that is configured to automatically read a probe identification code from a medical probe, including the generator 12, the generator quickly determines that the medical probe 302 has been previously used, and thus prevents the operation thereof (step 156).

The generator 12 periodically determines (e.g., every second, when operation of the medical probe 302 is based on the initial probe usage time T_i , and every time the medical probe 302, when operation of the medical probe 302 is based on the incremental probe usage U_i) whether the medical probe 302 has been previously used, and conditionally allows operation of the medical probe 302, while also conditionally disabling the probe identification component 50, based on this determination, until the medical probe 302 is disconnected from the generator 12 (steps 158-164). The control circuit 68 can periodically determine whether the medical probe 302 has been determined to have been previously used by periodically reading the probe identification code C_i from the probe identification

component 50, and comparing to the approved probe identification codes in the look-up list.

Alternatively, the generator 12 determines whether the medical probe 302 should be operated only one time per probe connection. That is, once operation of the medical probe 302 is allowed, the medical probe 302 can be operated without limitation until the medical probe 302 is physically disconnected from, and again connected to, the generator 12 (step 166).

Limiting Re-Use of Multiple Use Devices

Operation of the medical probe system 300 has been described above as limiting re-use of single use medical probes 302. The medical probe system 300, however, can be employed to limit re-usage of multiple use medical probes 302 as well. By way of non-limiting example, the medical probe system 300 can comprise a counter that keeps track of the number of times the medical probe 302 is indicated as being used. If the counter reaches a maximum limit, the medical probe 302 is prevented from being operated.

If the medical probe system 300 limits re-usage of the multiple use medical probe 302 based on the initial time that a medical probe 302 is therapeutically used, the previous initial probe usage time T_i is cleared and a current initial probe usage time T_i is stored each time the maximum predetermined elapsed time limit T_L expires. The counter, preferably located in the medical probe 302, keeps track of the number of times the maximum predetermined elapsed time limit T_L expires. The medical probe system 300 prevents further usage of the medical probe 302 once the maximum predetermined elapsed time limit T_L has expired and once the counter has reached a predetermined limit. For example,

if the number of uses of the medical probe 302 is to be limited to five, operation of the medical probe 302 may be prevented after the counter reaches five.

If the medical probe system 300 limits re-usage of a multiple use medical probe 302 based on detecting environmental changes, the probe sterilization indicator I_s is cleared

5 each time the medical probe 302 is therapeutically used and restored each time the medical probe 302 is sterilized. The counter, preferably located in the medical probe 302, keeps track of the number of times the presence of the probe sterilization indicator I_s is detected upon operation of the medical probe 302. The medical probe system 300 prevents further usage of the medical probe 302 once the presence of both the probe
10 sterilization indicator I_s and the previous probe usage indicator I_u is detected, and the counter has reached a predetermined limit.

Similarly, counters can be employed by the medical probe system 300 in limited re-usage of multiple use medical probes 302 based on cataloging probe usage, in which case multiple counters are preferably located in the RF generator 12 to respectively keep track
15 of a multiple number of medical probes 302. Or a counter can be employed by the medical probe system 300 in limiting re-usage of a medical probe 302 by preventing automatic identification of the medical probe 302, in which case the counter is preferably located in the medical probe 302.

It should further be noted that all of the above-described methods of limiting re-use
20 of a medical probe can each be employed alone to provide a single level safety measure, or alternatively, in combination, to provide a multiple level safety measure.

While preferred methods and embodiments have been shown and described, it will be apparent to one of ordinary skill in the art that numerous alterations may be made

without departing from the spirit or scope of the invention. Therefore, the invention is not to be limited except in accordance with the following claims.